

PRS-L70-F900-Si-PCB/CHP

Silicon piezo-resistive sensing cantilevers



General description


Piezo-Resistive Sensing (PRS) probes are silicon cantilevers with an integrated piezo-resistor for self-sensing scanning probe microscopy applications. The piezo-resistors are integrated into a matched Wheatstone bridge to raise the sensitivity and compensate environmental thermal drift. By using the self-sensing readout no laser adjustment is necessary in comparison to conventional optical readout AFM systems. This saves time during a cantilever change. The free space above the cantilever enables new applications and combination of AFM with various instruments. The cantilever chip is bonded to a small printed circuit board (PCB) with a small connector for a quick cantilever change. The counter part PCB for the cantilever PCB can be connected to a low-noise pre-amplifier with a flat flex cable.

Specifications	
Parameter	PRS-L70-F900-Si-PCB PRS-L70-F900-Si-CHP
Tip radius (apex)	<15 nm
Tip height	4...7 μm
Tip material	silicon
Resonant frequency	500...1300 kHz
Spring constant	35...400 N/m
AFM mode	tapping
sensitivity*	3 $\mu\text{V/nm}$
Length, width	70...85 μm , 30 \pm 1 μm
Material	silicon cantilever, boron doped 1k Ohm piezo resistors, aluminium tracks
Deflection sensing	on chip piezo-resistive bridge
Actuator	external shaker
Electrical connections	bonded to small PCB with connector (counter part PCB available) or optional bonding pads on chip
Chip dimensions (h, w, l)	0.3 / 1 / 2.8 mm
* not amplified, 2.048 V bridge supply	


Applications:

- Integration on a standard AFM scanner and high-speed AFM
- Force or deflection measurements within TEM, SEM, XPS, etc.


What about your application? Contact us!



Tip side of a PRS-L70 cantilever with Al tracks for reading out the sensor signal



Side view of a PRS-L70-F900 cantilever



Cantilever is bonded to a 6 x 4.5 mm PCB (height w. connector 1.6 mm, with CP-PCB: 2.5 mm); right: CP-flex-PCB



Hardware for amplified readout:
Low-noise pre-amplifier (45x35 mm)

SCL-Sensor.Tech. Fabrication GmbH

Viktor-Kaplan Straße 2, Bauteil E

2700 Wiener Neustadt, AUSTRIA

web: www.sclsensortech.com

Contact: Fabian Edlinger
Phone: +43 660 4424 871
fabian.edlinger@c-sense.at
Leaflet version: 2024-03-11

PRSA-L100-F500-Si-PCB

Silicon piezo-resistive sensing cantilevers



General description


Piezo-Resistive Sensing Active (PRSA) probes are silicon cantilevers with an integrated piezo-resistor and a heater for self-sensing and self-actuating scanning probe microscopy applications. The piezo-resistors are integrated into a matched Wheatstone bridge to raise the sensitivity and compensate environmental thermal drift. By using the self-sensing readout no laser adjustment is necessary in comparison to conventional optical readout AFM systems. This saves time during a cantilever change. The free space above the cantilever enables new applications and combination of AFM with various instruments. The cantilever chip is bonded to a small printed circuit board (PCB) with a small connector for a quick cantilever change. The counter part PCB for the cantilever PCB can be connected to a low-noise pre-amplifier with a flat flex cable.

Specifications	
Model	PRSA-L100-F500-Si-PCB
Tip radius (apex)	<15 nm
Tip height	4...6 µm
Tip material	silicon
Resonant frequency	200..800 kHz
Spring constant	8...530 N/m
Recomm. AFM mode	tapping, non-contact
sensitivity*	1...3 µV/nm
force sensitivity*	2.7...530 nN/µV
Length, width	110 ±5 µm, 48 ±2 µm
Material	silicon cantilever, boron doped 1k Ohm piezo resistors, aluminium tracks
Deflection sensing	on chip piezo-resistive bridge
Actuator	external shaker or on chip heater (12 +/-3 Ohm)
Electrical connections	bonded to small PCB with connector (counter part PCB available) or optional bonding pads on chip
Chip dimensions (h, w, l)	0.3 / 1.0 / 2.7 mm
* not amplified, 2.048 V bridge supply	


Applications:

- Integration on a standard AFM scanner and high-speed AFM
- Force or deflection measurements within TEM, SEM, XPS, etc.


What about your application? Contact us!




Tip side of a PRSA-L100 cantilever with Al tracks for reading out the deflection signal



Side view of a PRSA-L100 cantilever



Cantilever is bonded to a 6 x 4.5 mm PCB (height with connector 1.6 mm, with CP-PCB: 2.5 mm); left: counter part PCB



Hardware for amplified readout: Low-noise pre-amplifier (45x35 mm)

SCL-Sensor.Tech. Fabrication GmbH

Viktor-Kaplan Straße 2, Bauteil E
2700 Wiener Neustadt, AUSTRIA
web: www.sclsensortech.com

Contact: Fabian Edlinger
Phone: +43 660 4424 871
fabian.edlinger@c-sense.at
Leaflet version: 2024-03-11

PRSA-L300-Fxx-Si-PCB/CHP

Silicon piezo-resistive sensing cantilevers

xx...50/60/80



General description

Piezo-Resistive Sensing Active (PRSA) probes are silicon cantilevers with on chip integrated piezo-resistors and a heater for self-sensing and self-actuating scanning probe microscopy applications. The piezo-resistors are integrated into a matched Wheatstone bridge to raise the sensitivity and compensate environmental thermal drift. By using the self-sensing readout no laser adjustment is necessary in comparison to conventional optical readout AFM systems. This saves time during a cantilever change. The free space above the cantilever enables new applications and combination of AFM with various instruments. The cantilever chip is bonded to a small printed circuit board (PCB) with a small connector for a quick cantilever change. The counter part PCB for the cantilever PCB can be connected to a low-noise pre-amplifier with a flat flex cable.


Specifications

Model **	PRSA-L300-F60-Si-CHP	
	PRSA-L300-F50-Si-PCB	PRSA-L300-F80-Si-PCB
Tip radius (apex)	<15 nm	
Tip height	4...6 μm	
Tip material	silicon	
Resonant frequency **	30..65 kHz	65...95 kHz
Spring constant **	1...15 N/m	15...56 N/m
AFM mode	contact, non-contact	non-contact
sensitivity*	1...2 μV/nm	
force sensitivity*	0.5...56 nN/μV	
Length, width	300 ±5 μm, 110 ±3 μm	
Material	silicon cantilever, boron doped 1k Ohm piezo resistors, aluminium tracks	
Deflection sensing	on chip piezo-resistive bridge	
Actuator	external shaker or on chip heater (20-45 Ohm)	
Electrical connections	bonded to small PCB with connector (counter part PCB available) or optional bonding pads on chip	
Chip dimensions (h, w, l)	0.3 / 1.2 / 2.5 mm	
* not amplified (signal direct at the chip), 2.048 V bridge supply		
** Cantilever models are divided in two parameter ranges when electrical characterization is possible with bonded cantilevers		


Applications:

- Integration on a standard AFM scanner and high-speed AFM
- Force or deflection measurements within TEM, SEM, XPS, etc.


What about your application? Contact us!



Tip side of a PRSA-L300 cantilever with Al tracks for reading out the sensor signal



Side view of a PRSA-L300 cantilever



Cantilever is bonded to a 6 x 4.5 mm PCB (height with connector 1.6 mm, with CP-PCB: 2.5 mm); left: counter part PCB



Hardware for amplified readout:
Low-noise pre-amplifier (45x35 mm)

SCL-Sensor.Tech. Fabrication GmbH

Viktor-Kaplan Straße 2, Bauteil E

2700 Wiener Neustadt, AUSTRIA

web: www.sclsensortech.com

Contact: Fabian Edlinger
Phone: +43 660 4424 871
fabian.edlinger@c-sense.at
Leaflet version: 2024-03-11

PRSA-L400-F30-Si-PCB/CHP

Silicon piezo-resistive sensing cantilevers



General description


Piezo-Resistive Sensing Active (PRSA) probes are silicon cantilevers with integrated piezo-resistors on-chip and a heater for self-sensing and self-actuating scanning probe microscopy applications. The piezo-resistors are integrated into a matched Wheatstone bridge to raise the sensitivity and compensate environmental thermal drift. By using the self-sensing readout no laser adjustment is necessary in comparison to conventional optical readout AFM systems. This saves time during a cantilever change. The free space above the cantilever enables new applications and combination of AFM with various instruments. The cantilever chip is bonded to a small printed circuit board (PCB) with a small connector for a quick cantilever change. The cantilever PCB can be connected to a low-noise pre-amplifier by using our flexible counter-part PCB.

Specifications	
Model*	PRSA-L400-F30-Si-PCB PRSA-L400-F30-Si-CHP
Tip radius (apex)	<15 nm
Tip height	4.5...6.5 μm
Tip material	silicon
Resonant frequency	15...40 kHz
Spring constant (calc.)	0.4...8.3 N/m
Recomm. AFM mode	contact
sensitivity**	1 $\mu\text{V}/\text{nm}$
force sensitivity**	0.4...8.3 nN/ μV
Length, width	410 \pm 5 μm , 115 \pm 3 μm
Material	silicon cantilever, boron doped 1k Ohm piezo resistors, aluminium tracks
Deflection sensing	on chip piezo-resistive bridge
Actuator	external shaker or on-chip heater (11 +/-2 Ohm)
Electrical connections	bonded to small PCB with connector (counter part PCB available) or delivered as bare chip (bonding pads on-chip)
Chip dimensions (h, w, l)	0.3 / 1.2 / 2.5 mm
* Electrical characterization is only possible for bonded cantilevers.	
** measured without amplification, 2.048 V bridge supply	


Applications:

- Integration on a standard AFM scanner and high-speed AFM
- Force or deflection measurements within TEM, SEM, XPS, etc.


What about your application? Contact us!




Tip side of a PRSA-L400 cantilever with Al tracks for reading out the sensor signal



Side view of a PRSA-L400 cantilever



Cantilever is bonded to a 6 x 4.5 mm PCB (height with connector 1.6 mm, complete height connected to CP-PCB: 1.8 mm); left: counter part PCB



Hardware for amplified readout:
Low-noise pre-amplifier (45x35 mm)

SCL-Sensor.Tech. Fabrication GmbH

Viktor-Kaplan Straße 2, Bauteil E
2700 Wiener Neustadt, AUSTRIA
web: www.sclsensortech.com

Contact: Fabian Edlinger
Phone: +43 660 4424 871
fabian.edlinger@c-sense.at
Leaflet version: 2024-03-11